




AEGIS Liners

Sizing System



To determine proper BK liner sizing, two measurements are required:


- 1 Take circumference measurement 2" (5-6 cm) proximal to distal end of residual limb.
- 2 Take length measurement from mid-patella tendon to distal end of residual limb.

For liner size and item #, see chart below. Please specify Locking or Non-Locking (cushion) liner.

NOTE: On Aegis Streamline Liners, to specify the Covered Liner, add the letter "C" to the end of the item #.
All Aegis Ultimate, TF & SoFlex Liners come with cover.

AEGIS Ultimate				AEGIS Streamline				AEGIS SoFlex	
1 Limb Cir.	2 Length	Ultimate Locking	Ultimate Non-Locking	1 Limb Cir.	2 Length	Streamline Locking	Streamline Non-Locking	SoFlex Locking	SoFlex Non-Locking
14-20 cm	4"	BK U-XS-4L	BK U-XS-4C	17-19 cm	4"	BK 14-4S	BK 14-4G	SBK 14-4L	SBK 14-4C
	6"	BK U-XS-6L	BK U-XS-6C		6"	BK 14-6S	BK 14-6G	SBK 14-6L	SBK 14-6C
	8"	BK U-XS-8L	BK U-XS-8C		8"	BK 14-8S	BK 14-8G	SBK 14-8L	SBK 14-8C
20-26 cm	4"	BK U-S-4L	BK U-S-4C	19-21 cm	4"	BK 16-4S	BK 16-4G	SBK 16-4L	SBK 16-4C
	6"	BK U-S-6L	BK U-S-6C		6"	BK 16-6S	BK 16-6G	SBK 16-6L	SBK 16-6C
	8"	BK U-S-8L	BK U-S-8C		8"	BK 16-8S	BK 16-8G	SBK 16-8L	SBK 16-8C
26-32 cm	4"	BK U-M-4L	BK U-M-4C	21-23 cm	4"	BK 18-4S	BK 18-4G	SBK 18-4L	SBK 18-4C
	6"	BK U-M-6L	BK U-M-6C		6"	BK 18-6S	BK 18-6G	SBK 18-6L	SBK 18-6C
	8"	BK U-M-8L	BK U-M-8C		8"	BK 18-8S	BK 18-8G	SBK 18-8L	SBK 18-8C
32-38 cm	4"	BK U-L-4L	BK U-L-4C	23-25 cm	4"	BK 20-4S	BK 20-4G	SBK 20-4L	SBK 20-4C
	6"	BK U-L-6L	BK U-L-6C		6"	BK 20-6S	BK 20-6G	SBK 20-6L	SBK 20-6C
	8"	BK U-L-8L	BK U-L-8C		8"	BK 20-8S	BK 20-8G	SBK 20-8L	SBK 20-8C
38-42 cm	4"	BK U-XL-4L	BK U-XL-4C	10"	BK 20-10S	BK 20-10G	SBK 20-10L	SBK 20-10C	
	6"	BK U-XL-6L	BK U-XL-6C	25-27 cm	4"	BK 22-4S	BK 22-4G	SBK 22-4L	SBK 22-4C
	8"	BK U-XL-8L	BK U-XL-8C		6"	BK 22-6S	BK 22-6G	SBK 22-6L	SBK 22-6C
			8"		BK 22-8S	BK 22-8G	SBK 22-8L	SBK 22-8C	
				10"	BK 22-10S	BK 22-10G	SBK 22-10L	SBK 22-10C	
				27-29 cm	4"	BK 24-4S	BK 24-4G	SBK 24-4L	SBK 24-4C
					6"	BK 24-6S	BK 24-6G	SBK 24-6L	SBK 24-6C
					8"	BK 24-8S	BK 24-8G	SBK 24-8L	SBK 24-8C
				10"	BK 24-10S	BK 24-10G	SBK 24-10L	SBK 24-10C	
				30-32 cm	4"	BK 26-4S	BK 26-4G	SBK 26-4L	SBK 26-4C
					6"	BK 26-6S	BK 26-6G	SBK 26-6L	SBK 26-6C
					8"	BK 26-8S	BK 26-8G	SBK 26-8L	SBK 26-8C
				10"	BK 26-10S	BK 26-10G	SBK 26-10L	SBK 26-10C	
				32-34 cm	4"	BK 28-4S	BK 28-4G	SBK 28-4L	SBK 28-4C
					6"	BK 28-6S	BK 28-6G	SBK 28-6L	SBK 28-6C
					8"	BK 28-8S	BK 28-8G	SBK 28-8L	SBK 28-8C
				10"	BK 28-10S	BK 28-10G	SBK 28-10L	SBK 28-10C	
				35-37 cm	4"	BK 30-4S	BK 30-4G	SBK 30-4L	SBK 30-4C
					6"	BK 30-6S	BK 30-6G	SBK 30-6L	SBK 30-6C
					8"	BK 30-8S	BK 30-8G	SBK 30-8L	SBK 30-8C
				10"	BK 30-10S	BK 30-10G	SBK 30-10L	SBK 30-10C	

AEGIS TF			
1 Limb Cir.	Liner Size	AEGIS TF Locking	AEGIS TF Non-Locking
32-36 cm	30 cm	TF-30-L	TF-30-C
36-40 cm	34 cm	TF-34-L	TF-34-C
40-44 cm	38 cm	TF-38-L	TF-38-C
44-48 cm	42 cm	TF-42-L	TF-42-C
48-52 cm	46 cm	TF-46-L	TF-46-C



To determine proper AK liner sizing, one measurement is required:

- 1 Take circumference measurement 2" proximal to distal end of residual limb.

Engineered Silicone Products, LLC.

We're thinking on a higher level.

75 Mill Street, Newton, NJ 07860

888.932.7377 973.300.5120 fax: 973.300.5125

sales@wearesp.com www.wearesp.com

